Note: Document originally drafted in the English language.
注释：文件最初用英语起草。

Product Description

The Leo housing is designed for Varex Imaging rotating anode inserts having 71 mm (2.8 inches) or 80 mm (3 inches) diameter targets.

IEC Classification Class 1

Weight, Approximate:
Housing & Tube 17.4 kg (38.5 lbs)

Mounting Port plate or Trunnion
Metric Housing - M6-1 screws

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

Product Description	
Maximum Peak Voltage	150 kV
Cathode to Ground	75 kV
Anode to Ground	75 kV
Maximum X-Ray Tube Assembly Heat Content	900 kJ (1,250 kHU)
Nominal Continuous Input Power (maximum housing temperature 78°C)	200 W (278 HU/sec) IEC 60613:2010 @ sea Level
X-Ray tube assembly cooling is provided by ambient air flow.	
X-Ray Tube Assembly (Insert and Housing) Permanent Filtration ... 0.7mm Al/75kV IEC 60522/1999	
Loading Factors for Leakage Radiation 150 kV, 2.0 mA	
Temperature Limits for Storage and Transport	-20°C to +75°C
Humidity	10% to 90%
Atmospheric Pressure Range	70 kPa to 106 kPa
Thermal Switch	Normally Closed Contact Rating - 10A @ 240Vac
Open	80°C ±3.0°C (176°F ±6°F)
Federal Standard High Voltage Receptacles	(Complies to IEC 60526) NEMA Standard XR7-1979 (R1984, 1990)
X-Ray Tube Assembly (Complies to) IEC 60601-2-28	

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Dimensions are for Reference only
维度是供仅参考


Note 1

Focal Spot to Port Distance
80 mm Anode (RAD-14 & 74) [53.5 ±1.0] 2.11 ±0.040
71 mm Anode (RAD-12) [56.5 ±1.0] 2.23 ±0.040

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


Note: Window aperture is .827" (21mm) square
注释: 窗口开口是.827" (21mm)正方形

Cross Section of Port
端口的横截面

Stator Ratings and Characteristics
定子额定值及特征

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Stator Cord:

A six wire shield-grounded stator cord should be connected from the stator terminal to the motor control. All listed voltage and amperage measurements should be taken at the housing end.

If the stator cord is not provided, the cord used must be UL listed or a UL recognized component for all products used in the USA. Outside the USA, the cord used must comply with all applicable regulatory and statutory requirements for electrical and safety. The cord must have the following characteristics: six conductor, 18AWG (0.9mm²) Cu, rated 600V, 90°C minimum, shielded cable with PVC jacket 0.032 inch (0.081mm) thick minimum. Lead terminations are provided and must be used with wire that matches the above description.

Stator Power:

Time to full speed of the anode is a function of the power rating of the “starter” and the weight / diameter of the anode. All Varex Imaging stator types are rated for regular speed and high speed starters. Time to full speed for 71mm (2.8 inch) and 80mm (3.0 inch) anode series tubes is between 1.3 and 2.0 seconds.

Immediately following high speed anode rotation, the rotor speed must be reduced to 4000 r/min or less within 10 seconds using a suitable dynamic braking device.

No more than two high speed starts per minute are permissible. The starting voltage must never exceed 400 volts rms.

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
X-Ray Tube Assembly Heating and Cooling Curve
X 射线管汇编热化/冷却的曲线

Note:
Heat inputs into housing include tube power, filament power, and stator power. Heating curves based on no restrictions of natural convection around tube housing assembly.

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