Note: Document originally drafted in the English language.  
注释：文件最初用英语起草。

<b>Product Description</b>	
The Megalix-W™ housing is used for rotating anode inserts having 105mm (4 inch) and 120mm (4.75 inch) diameter targets.	
IEC Classification .....	Class 1
Weight, Approximate: Housing & Tube .....	26 kg (57 lbs)
Mounting .....	Trunnion
™All trademarks property of the respective manufacturer.	

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Stator - Wiring Diagram - Refer to the equipment manufacturer's instructions.

**Stator Power:**

Time to full speed of the anode is a function of the power rating of the "starter" and the weight / diameter of the anode. The Megalix stator is rated for regular speed and high speed starters. Time to full speed (2800 RPM) of 4.75 inch (120mm) anode series tubes is approximately 1 second.

Immediately following high speed anode rotation, the rotor speed must be reduced to 4000 r/min or less within 10 seconds using a suitable dynamic braking device.

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
X-Ray Tube Assembly Heating and Cooling Curve  
X 射线管汇编热化/冷却的曲线

**Note:**  
Heat inputs into housing include tube power, filament power, and stator power. Heating curves based on no restrictions of natural convection around tube housing assembly.

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