Note: Document originally drafted in the English language.

Product Description	
The Megalix-W™ housing is used for rotating anode inserts having 105mm (4 inch) and 120mm (4.75 inch) diameter targets.	
IEC Classification	Class 1
Weight, Approximate:	
Housing & Tube	26 kg (57 lbs)
Mounting	Trunnion
™All trademarks property of the respective manufacturer.	

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Product Description	
Maximum Potential Difference	125 kV
Cathode to Ground	62 kV
Anode to Ground	62 kV
Maximum X-Ray Tube Assembly Heat Content	1,800 kJ (2,400 kHU)
Maximum Continuous Heat Dissipation @25°C ambient (max. housing temperature 78°C) with Heat Exchanger	4,000 W (5,400 HU/sec)
X-Ray Tube Assembly (includes tube and housing) Permanent Filtration	1.5 mm Al @ 75 kV IEC 60522
Loading Factors for Leakage Radiation	125 kV, 4.0 mA
Temperature Limits for Storage and Transport	-20°C to +75°C
Humidity	10% to 90%
Atmospheric Pressure Range	70 kPa to 106 kPa
Pressure Switch	Normally Open 0.5 A @ 250 Vac
Close	8.5 PSIG (±1.5 PSIG)
Open	6.5 PSIG (±1.5 PSIG)
X-Ray Tube Assembly (Complies to)	IEC 60601-2-28

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Stator - Wiring Diagram - Refer to the equipment manufacturer's instructions.

Stator Power:

Time to full speed of the anode is a function of the power rating of the "starter" and the weight / diameter of the anode. The Megalix stator is rated for regular speed and high speed starters. Time to full speed (2800 RPM) of 4.75 inch (120mm) anode series tubes is approximately 1 second.

Immediately following high speed anode rotation, the rotor speed must be reduced to 4000 r/min or less within 10 seconds using a suitable dynamic braking device.

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