Note: Document originally drafted in the English language.

Product Description

The OPTI150™ housing is used for rotating anode inserts having 105mm (4 inch) and 120mm (4.75 inch) diameter targets.

IEC Classification. Class 1

Weight, Approximate:
Housing & Tube. 26 kg (57 lbs)

Mounting Trunnion

™All trademarks property of the respective manufacturer.

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Classification CEI Classe 1

Poids, Approximatif:
Gaine et Tube 26 kg (57 lbs)

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Product Description	
Maximum Potential Difference	150 kV
Cathode to Ground	75 kV
Anode to Ground	75 kV
Maximum X-Ray Tube Assembly Heat Content	1,800 kJ (2,400 kHU)
Maximum Continuous Heat Dissipation (max. housing temperature 85°C) without Heat Exchanger	300 W (400 HU/sec)
X-Ray Tube Assembly (includes tube and housing) Permanent Filtration	1.0 mm Al @ 75 kV IEC 60522
Loading Factors for Leakage Radiation	150 kV, 4.0 mA
Temperature Limits for Storage and Transport	-20°C to + 85°C
Humidity	10% to 90%
Atmospheric Pressure Range	70 kPa to 106 kPa
Pressure Switch	Normally Open 0.5 A @ 250 Vac
Close	8.5 PSIG (±1.5 PSIG)
Open	6.5 PSIG (±1.5 PSIG)
X-Ray Tube Assembly (Complies to)	IEC 60601-2-28
Tube assembly consists of housing, stator coil, x-ray tube (including anode, cathode) and insulating oil.	

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Stator Power:

Time to full speed of the anode is a function of the power rating of the "starter" and the weight / diameter of the anode. The stator is rated for regular speed and high speed starters. Time to full speed (2800 ±10% RPM) of 4.75 inch (120mm) anode series tubes is approximately 1 second.

Voltage: **Start** **Run**
 50/60Hz 250 VAC 75 VAC

Immediately following high speed anode rotation, the rotor speed must be reduced to 4000 r/min or less within 10 seconds using a suitable dynamic braking device.

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