Note: Document originally drafted in the English language.

**Product Description**

The B-135H housing is designed for Varex Imaging rotating Anode inserts having 102 mm (4 inch) diameter targets.

IEC Classification .....	Class 1
Weight, Approximate:	
Housing .....	20.4 kg (45 lbs)
With Trunnion Rings .....	22.9 kg (50.5 lbs)
With Trunnion Rings and Fan Assembly .....	24.1 kg (53.2 lbs)
Stator Cord .....	6.1 M (20 Feet)
Mounting .....	Trunnion or Port (using metric M6 screws)

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**Product Description**

Nominal X-ray Tube Voltage ..... 150 kV  
 Maximum Cathode to Ground ..... 75 kV  
 Maximum Anode to Ground ..... 75 kV  
 Grid to Cathode (If applicable) ..... -4 kV

Grid Control Voltages  
 Typical Bias Voltage for Cutoff at 150 kV ..... -3600 Vdc  
 Grid Voltage for Exposure ..... 0 Vdc

Maximum X-ray Tube Assembly Heat Content ..... 1480 kJ (2,000 kHU)

Nominal Continuous Input Power @ 25°C ambient  
 (max. housing temperature 78°C)  
 without air circulator ..... 222 Watts (300 HU/sec) IEC 60613:2010  
 with air circulator ..... 445 Watts (600 HU/sec) IEC 60613:2010

X-ray tube assembly cooling is provided by ambient air flow.

X-Ray Tube Assembly  
 Permanent Filtration ..... 0.7mm Al/75 kV IEC 60522/1999

High Voltage Cable Receptacles ..... Per IEC 60526

Loading Factors for Leakage Radiation ..... 150 kV, 4 mA

Temperature Limits for Storage and Transport ..... -20°C to +75°C  
 Humidity ..... +10% to +90%  
 Atmospheric Pressure Range ..... 70 kPa to 106 kPa

Thermal Switch ..... Normally Closed  
 7A @ 120Vac or 30 Vdc max.  
 Open ..... 85°C ±3.9°C (185°F ±7°F)

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


**Stator Power:**  
 Time to full speed of the anode is a function of the power rating of the "starter" and the weight / diameter of the anode. All Varex Imaging stator types are rated for regular speed and high speed starters. Time to full speed of 102 mm (4 inch) anode series tubes is between 1.3 and 2.0 seconds.  
 Immediately following high speed anode rotation, the rotor speed must be reduced to 4000 r/min or less within 10 seconds using a suitable dynamic braking device.  
 No more than two high speed starts per minute are permissible. The starting voltage must never exceed 600 volts rms.

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**Caution:**  
 Housings configured with heat exchangers having long hoses (> 5 ft or 1.5 m) are not configured for operation separate from the associated heat exchanger. In such cases, the tube unit must not be operated unless properly coupled to the heat exchanger. Reference the specific heat exchanger's product data sheet for the associated cooling curves.

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**Note:**  
 Heat inputs into housing include tube power, filament power, and stator power. Heating curves based on no restrictions of natural convection around tube housing assembly.

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