Note: Document originally drafted in the English language.

Product Description

The B-165H housing is designed for Varex Imaging rotating Anode inserts having 108 mm (4.25 inch) diameter targets.

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Soporte Rotable ó de entrada

Product Description

Nominal X-ray Tube Voltage	150 kV
Maximum Cathode to Ground	75 kV
Maximum Anode to Ground	75 kV
Grid to Cathode (If applicable)	-4 kV
Grid Control Voltages	
Typical Bias Voltage for Cutoff at 125 kV	-3700 Vdc
Grid Voltage for Exposure	0 Vdc
Maximum X-ray Tube Assembly Heat Content	1,100 kJ (1,500 kHU)
Nominal Continuous Input Power (max. housing temperature 78°C)	
without air circulator	200 Watts (270 HU/sec) IEC 60613:2010
with air circulator	400 Watts (540 HU/sec) IEC 60613:2010
with HE-100	1000 Watts (1400 HU/sec) IEC 60613:2010
X-ray tube assembly cooling is provided by ambient air flow.	
X-Ray Tube Assembly	
Permanent Filtration	1.0mm Al IEC 60522/1999
High Voltage Cable Receptacles	Per IEC 60526
Loading Factors for Leakage Radiation	150 kV, 4 mA
Temperature Limits for Storage and Transport	
Humidity	-20°C to +75°C
Atmospheric Pressure Range	+10% to +90%
	70 kPa to 106 kPa
Thermal Switch	
	Normally Closed
Open	3 A @ 120 Vac or 30 Vdc max. 76.7°C ±2.8°C (170°F ±5°F)
X-Ray Tube Assembly (Complies to)	IEC 60601-2-28

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
NOTE:

- REFERENCE DIMENSIONS ARE FOR INFORMATION PURPOSES AND ARE NOT REQUIRED FOR INSPECTION
- APERTURE SIZE OPTIONAL:
0.937 Sq. in. 1.250 Sq. in.
- CABLE RECEPTACLE
2 1/2 -20 UNF -2B THREAD

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
Stator Power:
 Time to full speed of the anode is a function of the power rating of the "starter" and the weight / diameter of the anode. All Varex Imaging stator types are rated for regular speed and high speed starters. Time to full speed of 102 mm (4 inch) anode series tubes is between 1.3 and 2.0 seconds.
 Immediately following high speed anode rotation, the rotor speed must be reduced to 4000 r/min or less within 10 seconds using a suitable dynamic braking device.
 No more than two high speed starts per minute are permissible. The starting voltage must never exceed 600 volts rms.

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