Rotating Anode X-Ray Tube Housing
 Boîtier de tube à rayons X à anode tournante
 Gehäuse der Röntgenröhre mit rotierender Anode
 Cubierta del Tubo de Rayos X de Ánodo Giratorio

Note: Document originally drafted in the English language.

Product Description

The MAM-RAD 100H housing is designed for Varex Imaging rotating anode mammography inserts having 100 mm (4 inch) diameter targets.

Housing: Shockproof, Aluminum; Leadlined; filled under vacuum with high dielectric insulating oil and hermetically sealed. The housing contains a stator for driving the rotor and fittings for using external Heat Exchanger.

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Tube Housing Specifications

Maximum Peak Voltage	
Anode to Cathode	49.9 kV
Anode to Ground	49.9 kV
Cathode to Ground	0 kV
Thermal Characteristics:	
Housing heat storage capacity	700,000 HU
Nominal Continuous Input Power	185 Watts
With HE50	600 Watts
With HE100	1130 Watts
IEC Classification Class 1	
Weight: 35 lbs. (16 kg) without cables or accessories	
X-Ray Tube Assembly	
Permanent Filtration	0.76 mm beryllium standard IEC 60522
(Optional molybdenum filters available) (0.0 mm Al/49.9 kV)	
Leakage Technique Factor	50 kV, 6.0 mA
Thermal Switch	
Normally Closed	Rating - 10A @ 240VAC
Open	60°C ±3°C (140°F ±5°F)
Closed	49°C ±3°C (120°F ±6°F)
Ambient Air Temperature Limits for Operation 5°C to 40°C	
Temperature Limits for Storage and Transport -20°C to +75°C	
Humidity	+10% to +90%
Atmospheric Pressure Range	70 kPa to 106 kPa

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
Stator: The MAM-RAD 100H housing is equipped with a standard Varex Imaging stator requiring a 24-30 microfarad phase shifting capacitor for 50/60 Hz., 2,850/3,450 rpm operation.

Stator Coil Resistance: Pin 1-3, 50 ohms; Pin 2-3, 20 ohms; Pin 1-2, 70 ohms.

Stator Motor Ratings:		Pin 3
Frequency (Hz.)	Volts	Amperes (Typical)
50/60Hz	50-60*	2.0
	120	4.2
	220	8.2

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