Note: Document originally drafted in the English language.

Product Description

The Sapphire housing is designed for Varex Imaging RAD series rotating anode inserts having 100mm (4 inch) diameter targets.

IEC Classification Class 1

Weight, Approximate:

Housing & Tube 18.9 kg (41.6 lbs)

Mounting Port plate or Trunnion
Standard Housing - 1/4-20 Screws
Metric Housing - M6-1 Screws

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Product Description

Maximum Potential Difference 150 kV
 Cathode to Ground 82.5 kV
 Anode to Ground 82.5 kV

Maximum X-Ray Tube Assembly Heat Content 1,111 kJ (1,500 kHU)

Maximum Continuous Heat Dissipation
 (max. housing temperature 78°C)
 without air circulator 216 W (300 HU/sec)
 with air circulator 432 W (600 HU/sec)

X-ray tube assembly cooling is provided by ambient air flow.

X-Ray Tube Assembly (Insert and Housing)
 Permanent Filtration 0.7mm Al/75kV IEC 60522

Loading Factors for Leakage Radiation 150 kV, 4.0 mA

Humidity +10% to +90%
 Temperature Limits for Storage and Transport -20°C to +75°C
 Atmospheric Pressure Range 70 kPa to 106 kPa

Thermal Switch Normally Closed
 contact rating - 7A @ 120Vac or 30 Vdc Max.
 Open 85°C ±3.9°C (185°F ±7°F)

Federal Standard High Voltage Cable Complies to IEC 60526

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


Filter Assembly
 Montage de Filtre
 Filteraubau
 Filtro Asamblado

Optional
 Filter Kit - IK163
 P/N 7576

Note: Window aperture is .875" (22mm) square on 16° mode and .656" (16.7mm) square on 10° and 12° models.

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Stator Power:

Time to full speed of the anode is a function of the power rating of the "starter" and the weight / diameter of the anode. All Varex Imaging stator types are rated for regular speed and high speed starters. Time to full speed of 102 mm (4 inch) anode series tubes is between 1.3 and 2.0 seconds.

Immediately following high speed anode rotation, the rotor speed must be reduced to 4000 r/min or less within 10 seconds using a suitable dynamic braking device.

No more than two high speed starts per minute are permissible. The starting voltage must never exceed 600 volts rms.

Stator Cord:

A six wire shield-grounded stator cord should be connected from the stator terminal to the motor control. All listed voltage and amperage measurements should be taken at the housing end.

If the stator cord is not provided, the cord used must be UL listed or a UL recognized component for all products used in the USA. Outside the USA, the cord used must comply with all applicable regulatory and statutory requirements for electrical and safety. The cord must have the following characteristics: six conductor, 18AWG (0.9mm²) Cu, rated 600V, 90°C minimum, shielded cable with PVC jacket 0.032 inch (0.081mm) thick minimum. Lead terminations are provided and must be used with wire that matches the above description.

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