Rotating Anode X-Ray Tube
Tubes Radiogénés à Anode Tournante
Röntgenröhre mit rotierender Anode
Tubos de Rayos-X con Ánodo Giratorio

■ Large - Black
Grand - Noir
Gross - Schwarz
Largo - Negro


■ Small - White
Petit - Blanc
Klein - Weiss
Pequeño - Blanco

○ Stand - By
Attente
Bereit Stehen
En Espera

□ Frame or Chassis
Masse
Chassis
Soporte o Chasis

○ X-Ray Tube
Tube Radiogène
Röntgenröhre
Tubo de Rayos X

○ Radiation Filter or Filtration
Filtre de rayonnement
Filterung
Filtración de Radiación

Note: Document originally drafted in the English language.

Product Description

The A-142 is a 2.8" (71mm) 150 kV, 222 kJ (300 kHU) maximum anode heat content, rotating anode insert. This insert is specifically designed for general radiographic and fluoro/spotfilm procedures. The insert features a 10° rhodium-tungsten molybdenum target and is available with the following nominal focal spots:

0.3 - 1.0
IEC 60336

Nominal Anode Input Power

Small - 9 kW IEC 60613
Large - 70 kW IEC 60613

For the equivalent anode input power of 80 Watts

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ANGIOGRAPHIC RATINGS

HOW TO USE ANGIOGRAPHIC CHARTS

General: Serial Radiography puts a severe demand on the x-ray tube due to the large number of exposures made in rapid succession. Intervals between exposures are fixed and so short that it is not possible for the anode track to cool to any extent during the exposure series. Therefore, the temperature of the anode track increases from exposure to exposure. The kW values used in the angiographic charts have been determined to prevent damage to the anode. The angiographic rating charts are usable to 100% anode heat storage. Exceeding 100% anode heat storage will cause anode track erosion with high risk of tube destruction.

Definition of Terms

Number of Exposures in Series: The number of exposures made in succession or the number of exposures made during one contrast injection.

Exposure Rate: The number of exposures made per second. For a series of exposures where the exposure rate changes, it must be assumed that all exposures will be made at the maximum rate. For example, if during a series 10 exposures will occur at one per second and 30 exposures at 4 per second, use the kW ratings in the 40 exposure column at 4 per second rate.

Exposure Time: Time in seconds of each exposure.

USING THE CHARTS:

Select Correct Chart:

50/60 or 150/180 Hz

0.3 or 1.0 Focal Spot

Note: 150/180 Hz rotor speed recommended for all angiography.

Determine the number of exposures in Series: With cut film angiography the number of exposures are known, however in Digital Angiography the number of exposures commonly are not known. When determining the number of exposures, assume worst case or past history.

Note: Most angiographic x-ray tubes fail from underestimating the number of exposures made in a series.

Determine kW of each exposure in Series: Referring to chart –find block under “Number of Exposures in Series” that is greater than or equal to expected number of exposures in Series. On left side directly opposite this block under “Exposure Rate per Second” column, select maximum rate per second that will be used for the exposure series. At the intersection of exposure rate and exposure time in seconds, find maximum kW allowed for each exposure.

kW = pkV x mA: The kW of the exposure can be any combination of mA and pkV allowed by the Radiographic and Filament Emission charts.

For Example: 80 pkV and 500 mA = 40 kW

Example: From chart A-142 150/180 Hz 3 Phase

1.0 Focal Spot, determine kW allowed with following known factors.

Maximum number of exposures40

Exposure time .050 second (50 milliseconds)

Maximum Exposure per second4

From chart find 40 exposure block. On left side directly opposite this block under “Exposure Rate per Second” column, select 4 exposures per second. Find .050 seconds at top of chart. At intersection of exposure rate line and exposure time, find 39.2 kW.

0.3 Focal Spot 10 Degrees $3\varnothing$ 50/60 Hz
0.3 Dimension Focale 30° 10 Degrés 50/60 Hz
0.3 Brennfleck 30° 10 Grad 50/60 Hz
0.3 De Marcas Focales 30° 10 Grados 50/60 Hz

Serial Load Ratings IEC 60613
Abaques de charges successives CEI 60613
Serienbetrieb-Belastungskurven IEC 60613
Ratio de carga en serie IEC 60613

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Note:

- (kW) of Exposure Equals mA x kV. For Example: 70 kV x 300 mA = 21 kW.
- Exposures less than .010 seconds will have a kW rating same as .010 seconds.

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