Rotating Anode X-Ray Tube Tubes  
Radiogènes à Anode Tournante  
Drehanoden - Röntgenröhre  
Tubos de Rayos-X con Ánodo Giratorio

Note: Document originally drafted in the English language.  
Note : Document à l'origine rédigé dans l'anglais.  
Anmerkung: Dokument ursprünglich gezeichnet in der englischen Sprache.  
Nota: Documento elaborado originalmente en la lengua inglesa.

### Product Description

The GS-2072 is a 5.25" (133 mm) 150 kV, 1,429 kJ (2.0 MHU) maximum anode heat content, rotating anode insert. This insert is specifically designed for CT Scanner. The insert features a 7° tungsten-rhenium facing on molybdenum with a graphite backed target and is available with the following nominal focal spots:

0.7 x 1.3  
IEC 60336

**Loading Factor for slit focal:**  
120 kV, 200 mA

**Maximum Anode Cooling Rate:**  
4,000 W (5,600 HU/sec)

**Maximum continuous anode heat dissipation:**  
3,000 W (4,200 HU/sec)

**Nominal CT Anode Input Power:**  
Large - 24 kW IEC 60613:2010

**Nominal CT Scan Power Index:**  
Large - 24 kW IEC 60613:2010

**Reference Angle:**  
Perpendicular to port face.

This insert is intended for use in a Varex Imaging B-200H housings.

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**Note:**  
 1. Limits are based on maximum track rating except for the following codes:  
 a - Limited by available heat storage.  
 b - Limited by window heating.  
 c - Limited by filament emission.  
 2. H.S. = Heat Storage  
 kV = Tube Voltage

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**Note:**  
 Rating charts reflect maximim tube performance. Tube operation is ultimately limited by system software.

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
Three Phase Emission (± .15A)  
 GS-2072 0.7 x 1.3

Maximum Peak Voltage	150 kV
Anode to Ground	75 kV
Cathode to Ground	75 kV
Maximum X-ray Tube Assembly Heat Content	3.1 MJ (4.3 MHU)
Maximum Housing Temperature	78°C
Nominal Continuous Input Power (Includes stator heat)	3.2 kW (4,480 HU/sec) IEC 60613:2010
Maximum Heat Exchanger Dissipation	5.0 kW (7,050 HU/sec)
Focal Point Position (Central Ray) Within 1mm (X, Y Direction from the center of radiation port.)	
X-Ray Tube Assembly Permanent filtration	1.0 mm Al IEC 60522
Loading Factors for Leakage Radiation	150 kV, 20 mA
High Voltage Cable Receptacles	Per IEC 60526
Ambient Air Temperature Limits for Operation	5°C to 40°C
Temperature Limits for Storage and Transport	-20°C to +75°C
Humidity	+10% to +90%
Atmospheric Pressure Range	70 kPa to 106 kPa
Weight - Housing	31.5 kg (69.4 lbs)
Heat Exchanger	16.3 kg (35.9 lbs)
IEC Classification	Class 1
Safety Devices: Thermal Switch Normally closed contact	Opening at 85°C ±3.9°C
Pressure Switch Normally Open contact	closes before 5.0 PSID Rising opens at 3.0 PSID ±0.2 PSID Falling
Filament Frequency Limits	50 HZ - 25 KHZ
Power Supply	DC

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

## X-ray Tube Assembly Heating and Cooling Curve

**Note:**

- Heat input into housing includes all power sources; tube, filament, stator and circulating pump.
- Heating curves based on no restrictions to air flow through heat exchanger, or natural convection around tube housing assembly.
- Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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

<b>Stator Type:</b> "R" Stator		
<b>Stator Coil Resistance:</b>		
Black to White	14.0 Ohms ±15%	
Green to White	46.0 Ohms ±15%	
<b>Starter Voltage:</b>		
50/60 Hz	Start 200 VAC	Run 70 VAC
<b>Time to Full Speed:</b>		
50/60 Hz	0 - 2700 RPM	17 Sec.
<b>X-Ray Tube Assembly:</b>		
GS-2072/B-200H	IEC 60601-2-28	

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**Note:**  
 1. Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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