Note: Document originally drafted in the English language.
注释: 文件最初用英语起草。

Product Description

The GS-3070 is a 5.5" (140 mm) 150 kV, 2.5 MJ (3.5 MHU) maximum anode heat content, rotating anode insert. This insert is specifically designed for CT Scanners. The insert features a 7° tungsten-rhenium facing on molybdenum with a graphite backed target and is available with the following nominal focal spot:

0.7 x 0.8
1.2 x 1.4
IEC 60336

Loading Factor for slit focal:

Small - 120 kV, 100 mA
Large - 120 kV, 200 mA

Maximum Anode Cooling Rate:

8,750 W (12,250 HU/sec)

Continuous Anode Input Power:

3,400 W (4,760 HU/sec)

Nominal CT Anode Input Power:

Small - 21 kW IEC 60613
Large - 45 kW IEC 60613

Nominal CT Scan Power Index:

Small - 21 kW IEC 60613:2010
Large - 41 kW IEC 60613:2010

Reference Axis:

Perpendicular to port face.

This insert is intended for use in a Varex Imaging B-240H housing.

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

Note:
 1. Limits are based on maximum track rating except for the following codes:
 a - Limited by available heat storage.
 b - Limited by window heating.
 c - Limited by filament emission.
 2. H.S. = Heat Storage
 kV = Tube Voltage

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Note:
 Rating charts reflect maximum tube performance. Tube operation is ultimately limited by system software.

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3 Ø ≡

Product Description






Maximum Peak Voltage	150 kV
Anode to Ground	75 kV
Cathode to Ground	75 kV
Maximum X-ray Tube Assembly Heat Content	3.8 MJ (5.2 MHU)
Nominal Continuous Input Power (max. housing temperature 78°C) (Includes stator heat)	3.7 kW (5.18 kHU/sec)
Maximum Heat Exchanger Dissipation	5.0 kW (7.0 kHU/sec)
Focal Point Position (Central Ray) Within 1mm (X, Y Direction from the center of radiation port.)	
X-Ray Tube Assembly	
Permanent filtration	1.0 mm Al IEC 60522
Loading Factors for Leakage Radiation	150 kV, 23 mA
High Voltage Cable Receptacles	Per IEC 60526
Ambient Air Temperature Limits for Operation	5°C to 40°C
Temperature Limits for Storage and Transport	-20°C to +75°C
Humidity	+10% to +90%
Atmospheric Pressure Range	70 kPa to 106 kPa
Weight: Housing	33.1 kg (73 lbs)
Heat Exchanger	16.3 kg (35.9 lbs)
IEC Classification	Class 1
Safety Devices: (Internal) Thermal Switch	
Normally Closed Contact	Opening at 85°C ±3.9°C
Pressure Switch	
Normally Open contact	closes before 5.0 PSID Rising opens at 3.0 PSID ±0.2 PSID Falling
Filament Frequency Limits	50 Hz - 25 kHz
Power Supply	DC

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
Dimensions are for reference only
维度是供仅参考

Dimensions are in millimeters
尺寸单位为毫米

Tube Housing Assembly Heating and Cooling
管套组件加热与冷却

Note:


- Heat input into housing includes all power sources; tube, filament, stator and circulating pump.
- Heating curves based on no restrictions to air flow through heat exchanger, or natural convection around tube housing assembly.
- Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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
Stator Type: "R"		
Stator Coil Resistance:		
Black to White	14 Ω	±15%
Green to White	46 Ω	±15%
Starter Voltage:		
	Start	Run
50/60 Hz	200 VAC	80 VAC
150/180 Hz	470 VAC	140 VAC
Time to Full Speed:		
50/60 Hz	0 - 3000 RPM	9.6 Sec.
150/180 Hz	0 - 8000 RPM	9.6 Sec.
X-Ray Tube Assembly:		
GS-3070/B-240H	IEC 60601-2-28	

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Note:
Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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