Rotating Anode X-Ray Tube  
 Tubes Radiogènes à Anode Tournante  
 Röntgenröhre mit rotierender Anode  
 Tubos de Rayos-X con Ánodo Giratorio

Note: Document originally drafted in the English language.

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**Note:**  
 1. Limits are based on maximum track rating except for the following codes:  
 a - Limited by available heat storage.  
 b - Limited by window heating.  
 c - Limited by filament emission.  
 2. H.S. = Heat Storage  
 kV = Tube Voltage

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**Note:**  
 Rating charts reflect maximum tube performance. Tube operation is ultimately limited by system software.

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
Le Gaine B-245H

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Encaje de B-245H

Maximum Peak Voltage	150 kV
Anode to Ground	75 kV
Cathode to Ground	75 kV
Maximum X-ray Tube Assembly Heat Content	1.5 MJ (2.0 MHU)
Maximum Continuous Heat Dissipation (max. housing temperature 78°C) (Includes stator heat)	3.25 kW (4.55 kHU/sec)
Maximum Heat Exchanger Dissipation	3.25 kW (4.55 kHU/sec)
Focal Point Position (Central Ray) Within 1mm (X, Y Direction from the center of radiation port.)	
X-Ray Tube Assembly	
Permanent filtration	1.0 mm Al IEC 60522
Loading Factors for Leakage Radiation	150 kV, 23 mA
Federal Standard High Voltage Cable	72
Ambient Air Temperature Limits for Operation	5°C to 40°C
Temperature Limits for Storage and Transport	-20°C to +75°C
Humidity	+10% to +90%
Atmospheric Pressure Range	70 kPa to 106 kPa
Weight: Housing	35 kg (77.1 lbs)
Heat Exchanger	17.3 kg (38.1 lbs)
IEC Classification	Class 1
Safety Devices: Thermal Switch	
Normally Closed Contact	Opening at 85°C ±3°C
Flow Switch - Normally Open contact	Contacts close with adequate oil flow.
Filament Frequency Limits	50 Hz - 25 kHz
Power Supply	DC

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

B-245H Housing  
 Le Gaine B-245H  
 Das B-245H Gehäuse  
 Encaje de B-245H

**Note:**

- Heat inputs into housing include tube power, filament power, and stator power.
- Heating curves based on no restrictions of natural convection around tube housing assembly.
- Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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

<b>Stator Type:</b> "R"
<b>Stator Coil Resistance:</b> Black to White 14 Ω ±15% Green to White 46 Ω ±15%
<b>Starter Voltage:</b> <u>Start</u> <u>Run</u> 50/60 Hz 230 VAC 85 VAC
<b>Time to Full Speed:</b> 50/60 Hz 0 - 2800 RPM 10 Sec.
<b>X-Ray Tube Assembly:</b> GS-3074/B-245H IEC 60601-2-28

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**Note:**  
 Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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