Rotating Anode X-Ray Tube Tubes
Radiogènes à Anode Tournante
Drehanoden - Röntgenröhre
Tubos de Rayos-X con Ánodo Giratorio

Note: Document originally drafted in the English language.

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Product Description

The GS-3075A is a 5.5" (140 mm) 150 kV, 2.5 MJ (3.5 MHU) maximum anode heat content, rotating anode insert. This insert is specifically designed for Shimadzu CT Scanners capable of 0.75 second one rotation scanning. The insert features a 7° tungsten-rhenium facing on molybdenum with a graphite backed target and is available with the following nominal focal spot:

0.7 x 1.2
1.2 x 1.5
IEC 60336

Loading Factor for slit focal:

Small - 120 kV, 200 mA
Large - 120 kV, 300 mA

Maximum Anode Cooling Rate:

8,750 W (12,250 HU/sec)

Maximum continuous anode heat dissipation:

3,400 W (4,750 HU/sec)

Nominal Anode Input Power:

Small - 25 kW IEC 60613
Large - 40 kW IEC 60613

Reference Axis:

Perpendicular to port face.

This insert is intended for use in Varex Imaging B-240H housing.

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Note:
 1. Limits are based on maximum track rating except for the following codes:
 a - Limited by available heat storage.
 b - Limited by window heating.
 c - Limited by filament emission.
 2. H.S. = Heat Storage
 kV = Tube Voltage

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Note:
 Rating charts reflect maximum tube performance. Tube operation is ultimately limited by system software.

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

Le Gaine B-240H

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Encaje de B-240H

Maximum Peak Voltage	150 kV
Anode to Ground	75 kV
Cathode to Ground	75 kV
Maximum X-ray Tube Assembly Heat Content	1.5 MJ (2.0 MHU)
Maximum Housing Temperature	78°C
Maximum Continuous Heat Dissipation (Includes stator heat)	3.7 kW (5.18 kHU/sec)
Maximum Heat Exchanger Dissipation	5.0 kW (7.0 kHU/sec)
Focal Point Position (Central Ray) Within 1mm (X, Y Direction from the center of radiation port.)	
X-Ray Tube Assembly	
Permanent filtration	1.0 mm Al IEC 60522
Loading Factors for Leakage Radiation	150 kV, 23 mA
Federal Standard High Voltage Cable	72
Ambient Air Temperature Limits for Operation	5°C to 40°C
Temperature Limits for Storage and Transport	-20°C to +75°C
Humidity	+10% to +90%
Atmospheric Pressure Range	70 kPa to 106 kPa
Weight: Housing	33.1 kg (73 lbs)
Heat Exchanger	15.9 kg (35 lbs)
IEC Classification	Class 1
Safety Devices: (Internal) Thermal Switch	
Normally Closed Contact	Opening at 85°C
Flow Switch - Normally Open contact	
	Contacts close with adequate oil flow.
Filament Frequency Limits	50 HZ - 25 kHz
Power Supply	DC

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

Note:
 1. Heat inputs into housing include tube power, filament power, and stator power.
 2. Heating curves based on no restrictions of natural convection around tube housing assembly.
 3. Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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

Stator Type: "R"	
Stator Coil Resistance: Black to White 14.0 Ohms ±15% Green to White 46.0 Ohms ±15%	
Starter Voltage:	Start 370 VAC Run 80 VAC
Time to Full Speed: 50/60 Hz	0 - 2700 RPM 6 Sec.
X-Ray Tube Assembly: GS-3075A/B-240H IEC 60601-2-28	

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Note:
 1. Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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