Note: Document originally drafted in the English language.
 Note : Document à l'origine rédigé dans l'anglais.
 Anmerkung: Dokument ursprünglich gezeichnet in der englischen Sprache.
 Nota: Documento elaborado originalmente en la lengua inglesa.

Product Description

The MCS-6200 is a 7.8" (200 mm) 140 kV, 4.3 MJ (6.0 MHU) maximum anode heat content, rotating anode insert. This insert is specifically designed for CT Scanners. The insert features a 7° tungsten-rhenium facing on molybdenum with a graphite backed target and is available with the following nominal focal spots:

0.7 x 0.8
 1.2 x 1.4
 IEC 60336

Loading Factor:

Small - 120 kV, 200 mA
 Large - 120 kV, 200 mA

Maximum Anode Cooling Rate:
 8,800 W (12,320 HU/sec)

Maximum continuous
 anode heat dissipation:
 8,800 W (12,320 HU/sec)

Nominal CT Anode Input Power:
 Small - 33 kW IEC 60613:2010
 Large - 60 kW IEC 60613:2010

Nominal CT Scan Power Index:
 Small - 30 kW IEC 60613:2010
 Large - 59 kW IEC 60613:2010

Reference Axis:
 Perpendicular to port face.

This insert is intended for use in
 Varex Imaging B-61XH housing.

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Note:
 1. Limits are based on maximum track rating except for the following codes:
 a - Limited by available heat storage.
 b - Limited by window heating.
 c - Limited by filament emission.
 2. H.S. = Heat Storage
 kV = Tube Voltage

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
Note:
 Rating charts reflect maximim tube performance. Tube operation is ultimately limited by system software.

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
THREE PHASE EMISSION (± .15 A)
 0.7 x 0.8 

THREE PHASE EMISSION (± .15 A)
 1.2 x 1.4 



Maximum Peak Voltage	140 kV
Anode to Ground	0 kV
Cathode to Ground	140 kV
Maximum X-ray Tube Assembly Heat Content	7.0 MJ (9.8 MHU)
Nominal Continuous Input Power	4.0 kW (5.6 kHU/sec)
	IEC 60613:2010
Maximum Housing Temperature	78°C
Maximum Heat Exchanger Dissipation	6.0 kW (8.4 kHU/sec)
Focal Point Position (Central Ray) Within 1mm (X, Y Direction from the center of radiation port.)	
X-Ray Tube Assembly	
Permanent filtration	1.0 mm Al IEC 60522
Loading Factors for Leakage Radiation	140 kV, 29 mA
High Voltage Cable	Special
Ambient Air Temperature Limits for Operation	5°C to 40°C
Temperature Limits for Storage and Transport	-20°C to +75°C
Humidity	+10% to +90%
Atmospheric Pressure Range	70 kPa to 106 kPa
Weight - Housing	61.8 kg (136.3 lbs)
Heat Exchanger	32.5 kg (71.6 lbs)
IEC Classification	Class I
Safety Devices	
Housing	
Thermal Switch: Normally Closed Contacts	Opens at 85°C
Thermal Fuse: Normally Closed Contacts	Opens at 105°C
Heat Exchanger	
Flow Switch: Normally Open Contacts, closes with adequate oil flow	
Filament Frequency Limits	50 HZ - 40 KHZ
Power Supply	DC

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




Pin-outs for Low Voltage Feed-through

- A: Stator Phase Shift
- Pin B: Stator Main
- Pin C: Stator Common
- Pin D: Housing Ground
- E: Not Connected/Used
- Pins F & G: Thermal Switch

Low Voltage Feedthrough

TUBE HOUSING ASSEMBLY HEATING AND COOLING

Note:

- Heat input into housing includes all power sources; tube, filament, stator and circulating pump.
- Heating curves based on no restrictions to air flow through heat exchanger, or natural convection around tube housing assembly.
- Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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

Stator Type: "3 Ø"		
Stator Coil Resistance: 2.15 Ohms ± 15%		
Starter Voltage:	Start	Run
50 Hz	180-200 VAC	80 VAC
100 Hz	180-200 VAC	100 VAC
Time to Full Speed:		
50 Hz	10 Sec.	
100 Hz	10 Sec.	
X-Ray Tube Assembly: MCS-6200/B-61XH IEC 60601-2-28		

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ANODE HEATING AND COOLING CURVES

Note:

1. Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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