Rotating Anode X-Ray Tube Tubes
Radiogénés à Anode Tournante
Drehanoden - Röntgenröhre
Tubos de Rayos-X con Ánodo Giratorio

Note: Document originally drafted in the English language.
Note : Document à l'origine rédigé dans l'anglais.
Anmerkung: Dokument ursprünglich gezeichnet in der englischen Sprache.
Nota: Documento elaborado originalmente en la lengua inglesa.

Product Description

The MCS-8064 is a 9.4" (240 mm) 140 kV, 5.7 MJ (8.0 MHU) maximum anode heat content, rotating anode insert. This insert is specifically designed for GE CT Scanners. The insert features a 7° tungsten-rhenium facing on molybdenum with a graphite backed target and is available with the following nominal focal spots:

0.9 x 0.7
1.2 x 1.2
IEC 60336

Maximum Anode Cooling Rate:
21,800 W (30,520 HU/sec)

Maximum continuous anode heat dissipation:
8,000 W (11,200 HU/sec)

Nominal Anode Input Power:
Load Time 1 Second
Small - 47 kW IEC 60613:2010
Large - 100 kW IEC 60613:2010

Nominal CT Anode Input Power:
Small - 47 kW IEC 60613:2010
Large - 96 kW IEC 60613:2010

Nominal CT Scan Power Index:
Small - 40 kW IEC 60613:2010
Large - 81 kW IEC 60613:2010

Reference Axis:
Perpendicular to port face.

This insert is intended for use in a Varex Imaging housing B-680H.

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NOTE: 72 kW base configuration is not available in all markets.

Large Spot Single Exposure Limits				
Scan Time	140 kV	120 kV	100 kV	80 kV
5	615*	800*	760*	675†
10	575†	790*	760*	675†
20	535†	645†	760*	675†
30	490	560	695†	675†
40	445	510	635†	675†
50	410	485	585	675†
60	385	460	550	675†

* Available only with VCT 85 kW Option installed. Otherwise constrained to the max mA available.

* Available only with VCT 85 kW & VCT Hi Power Options installed. Otherwise constrained to the max mA available.

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The serial exposures are controlled by system software and are applicable for repeat every 10 minutes based on 3 hour wait period after tube warm up has completed.

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† Available only with VCT 85 kW Option installed. Otherwise constrained to the maximum mA available.

* Available only with VCT 85 kW & VCT Hi Power Options installed. Otherwise constrained to the maximum mA available.

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Highest Constant Load at 4s (Reference IEC 60601-2-44 6.8.2)

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The system can acquire 72 kW at 120 kVp for 4 seconds scan duration. The system can acquire 84 kW at 120 kVp for 4 seconds scan duration if the 85 kW Option is enabled. The system can acquire 96 kW at 120 kVp for 4 seconds scan duration if the 85 kW and Hi Power Options are enabled. The single exposures are controlled by system software.

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






MCS-8064

Note:

- Heat input into housing includes all power sources; tube, filament, stator and circulating pump.
- Heating curves based on no restrictions to air flow through heat exchanger, or natural convection around tube housing assembly.
- Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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