
Rotating Anode X-Ray Tube  
 Tubes Radiogènes à Anode Tournante  
 Röntgenröhre mit rotierender Anode  
 Tubos de Rayos-X con Ánodo Giratorio

- Common - Red  
 Neutre - Rouge  
 Neutral - Rot  
 Común - Rojo
- Large - Black  
 Grand - Noir  
 Gross - Schwarz  
 Largo - Negro
- Small -White  
 Petit - Blanc  
 Klein - Weiss  
 Pequeño - Blanco
- Stand - By  
 Attente  
 Bereitschaft  
 En Espera
- Frame or Chasis  
 Masse  
 Chassis  
 Soporte o Chasis
- X-Ray Tube  
 Tube Radiogène  
 Röntgenröhre  
 Tubo de Rayos X
- Radiation Filter or Filtration  
 Filtre de rayonnement  
 Filterung  
 Filtración de Radiación

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## CINERADIOGRAPHIC RATINGS

### HOW TO USE CINERADIOGRAPHIC CHARTS

**General:** With the Cineradiographic rating chart we can determine the maximum allowable kW of the Cine pulse, or with a given kW determine maximum time in seconds the Cine run can progress.

The Most common way of using the charts is to determine maximum time of any expected Cine run and maximum duty factor. With a known duty factor and Cine run time kW can easily be determined.

**Definition of Terms**

**Time in seconds:** Total time of one Cine run, usually 5 to 12 seconds.

**Duty Factor in Percent (DF%):** Actual time during one second the x-ray tube is producing x-rays. If we select a 4 msec pulse width and 60 exposures per second the x-ray tube will be producing x-rays for a total of 240 msec each second or 24% of the time. The higher the DF number, the more load placed on the x-ray tube.

**Peak Pulse Power:** Peak energy in watts of any one Cine Pulse. Can be any combination of kV and mA allowed by Radiographic and Filament Emission curves.

Example:           80 kV at 400 mA equals  
  
                          80,000 V x 0.4 A = 32,000 W or 32 kW

**USING THE CINE RATING CHARTS:**

RAD-94 150/180 Hz 3 Phase 0.8 Focal Spot

**Example:** Determine maximum kW allowed with the following known factors:  
Maximum Pulse Width ..... 4 msec  
Exposures per Second ..... 60  
Maximum Cine Run Time .....10 seconds

**Calculate Duty Factor:** (DF%)

$$DF\% = \frac{\text{Pulse Width (mSec)} \times \text{Frames per Second}}{10}$$

$$DF\% = \frac{4 \text{ msec} \times 60 \text{ exp/sec}}{10} = \frac{240}{10} = 24\%$$

Refer to Rating Chart RAD-94 150/180 Hz  
3 Phase 0.8 Focal Spot:

At bottom of chart find 10 second line. Move vertically to intersection with 24% DF curve. Make a horizontal reference to left side of rating chart and note kW rating of 34 kW.



We now know each pulse during the cine run can have a maximum rating of 34 kW under conditions given in example.

kW = kV x mA. The kW of the exposure can be any combination of mA and kV allowed by the Radiographic and Filament Emission Charts.



The Cine rating charts are usable to 100% anode heat storage. The start of Cine run should be below 70% and heat storage. Exceeding 100% anode heat storage will cause anode track erosion with high risk of tube destruction.

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

THREE PHASE EMISSION (± .15 A)


0.8 

Note:  
 When using these emission curves for trial exposures, refer to the power rating curves shown for maximum kV, tube emission, filament current, exposure time, and target speed.

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Salt Lake City, UT 1-801-972-5000

For a complete listing of our global offices,  
visit [www.vareximaging.com](http://www.vareximaging.com)

Manufactured by Varex Imaging Corporation  
Fabrique par Varex Imaging Corporation  
Hergestellt von Varex Imaging Corporation  
Fabricado por Varex Imaging Corporation

Specifications subject to change without notice.  
Spécifications susceptibles d'être modifiées sans préavis.  
Technische Daten ohne Gewähr.  
Especificaciones sujetas a cambio sin previo aviso.