

**DESCRIPTION:**

Heat Exchanger (oil to air) with single fan. For use with the B-330H housing. Oil pump is gear type.

**PERFORMANCE:**

2.4 kW with 50° delta T (194 kHU/minute). Refer to cooling chart for other delta T ratings.

**TOTAL HOSE LENGTH:**

30 inches minimum (0.75 meters)  
164 feet maximum (50 meters)

**POWER REQUIREMENTS:**

HE-321A 220 - 230 VAC 50/60 Hz 1.2 Amps

**FAIL SAFE DEVICE:**

Flow Switch - Normally open contacts.  
Flow Switch Operating Pressure - 1000PSIG  
Flow Switch Contact Rating - 250 Volts Max

**DIMENSIONS:**

Width - 13.75 inches (349 mm)  
Length - 16 inches (406 mm)  
Height - 5.25 inches (133 mm)

**WEIGHT:**

40 Pounds (18.1 kg)



**MOUNTING:**

Select mounting location where ambient temperature is low.

Connect power cord to source that is energized whenever x-ray tube with heat exchanger is selected.


Connect flow switch and thermal switch on housing in series with existing exposure interlocks.

**CONNECTING OIL LINES WITH QUICK DISCONNECT FITTINGS:**

To prevent air in housing from oil line quick disconnect fittings observe the following procedures.



- A. Do not apply any side pressure on fittings such as from tie-wrapping to cable bundles or looping in tight circle.
- B. Never run heat exchanger with fittings disconnected.
- C. **Always** return heat exchanger with x-ray tube regardless of reason for return.

**COOLING CHARACTERISTICS  
HE -321A POWER DISSIPATION**



### Outline Diagram

Dimensions are for reference only



DIMENSIONAL DATA		
	INCHES	mm
DIM	REF	REF
A	15.91	404.1
B	13.63	346.2
C	0.50	12.7
D	2.25	57.2
E	3.58	90.9
F	0.95	24.1
G	4.95	125.7
H	1.38	35.1
I	0.52	13.2
J	2.30	58.4
K	1.08	27.4
L	0.58	14.7
M	7.94	201.7
N	5.13	130.3
P	0.45	11.4
Q	10.13	257.3
R	6.27	159.3
S	10.88	276.4
T	5.35	135.9
U	4.19	106.4
V	1.33	33.8
W		
X		
Y		

4 X 1/4-20

