Note: Document originally drafted in the English language.

Product Description

The B-100 housing is designed for Varex Imaging rotating Anode inserts having 71 mm (2.8 inches) and 80 mm (3.0 inches) diameter targets.

IEC Classification Class 1
 Weight, Approximate:
 Housing 17.5 kg (38.5 lbs)
 Stator Cord 0.9m (3 Ft)
 Mounting Trunnion or Port Plate

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Product Description

Nominal X-ray Tube Voltage 150 kV
 Maximum Cathode to Ground 75 kV
 Maximum Anode to Ground 75 kV
 Grid to Cathode (if applicable) -4 kV

Grid Control Voltages
 Typical Bias Voltage for Cutoff at 150 kV -3600 Vdc
 Grid Voltage for Exposure 0 Vdc

Housing Heat Storage Capacity 926 kJ (1,250 kHU)

Nominal Continuous Input Power
 without air circulator 185 Watts (250 HU/sec) IEC 60613:2010
 with air circulator 370 Watts (500 HU/sec) IEC 60613:2010

X-Ray Tube Assembly
 Permanent Filtration 0.7 mm Al/75 kV IEC 60522/1999

Leakage Technique Factors 150 kV, 3.3 mA

Temperature Limits for Storage and Transport -20°C to +75°C
 Humidity +10% to +90%
 Atmospheric Pressure Range 70 kPa to 106 kPa

Thermal Switch 7A @ 120 VAC or 30 VDC max.
 Normally Closed Opens @ 85°C ±3.9°C (185°F ±7°F)

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




Receptacle Key
 Clef du Receptacle
 Hochspannungsbuchsen
 Llava del Receptaculo

Note
 Focal Spot to Port Distance
 80 mm Anode 2.1 inches (53 mm)
 71 mm Anode 2.2 inches (56 mm)

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
To ensure proper target rotation, each stator cord is equipped with a special twist-locking connector and a mating connector.

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

Stator Power:
 Time to full speed of the anode is a function of the power rating of the "starter" and the weight / diameter of the anode. All Varex Imaging stator types are rated for regular speed and high speed starters. Time to full speed for 71mm (2.8 inch) and 80mm (3.0 inch) anode series tubes is between 1.3 and 2.0 seconds.
 Immediately following high speed anode rotation, the rotor speed must be reduced to 4000 r/min or less within 10 seconds using a suitable dynamic braking device.
 No more than two high speed starts per minute are permissible. The starting voltage must never exceed 500 volts rms.

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Note:
Heat inputs into housing include tube power, filament power, and stator power. Heating curves based on no restrictions of natural convection around tube housing assembly.

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Salt Lake City, UT 1-801-972-5000

For a complete listing of our global offices,
visit www.vareximaging.com

Manufactured by Varex Imaging Corporation
Fabrique par Varex Imaging Corporation
Hergestellt von Varex Imaging Corporation
Fabricado por Varex Imaging Corporation

Specifications subject to change without notice.
Spécifications susceptibles d'être modifiées sans préavis.
Technische Daten ohne Gewähr.
Especificaciones sujetas a cambio sin previo aviso.