Note: Document originally drafted in the English language.

**Product Description**

The B-145A housing is designed for Varex Imaging rotating Anode inserts having 71 mm (3 inch) diameter targets.

IEC Classification ..... Class 1

Weight, Approximate:  
Housing & Tube ..... 14.51 kg (32 lbs)

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Classification CEI ..... Classe 1

Poids, Approximatif:  
Gaine et Tube ..... 14,51 kg (32 lbs)

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**Product Description**

Maximum Potential Difference ..... 125 kV  
 Cathode to Ground ..... 63 kV  
 Anode to Ground ..... 63 kV

Maximum X-Ray Tube Assembly Heat Content ..... 500 kJ (675 kHU)

Nominal Continuous Input Power  
 (Maximum Housing Temperature 78°C) ..... 140 W (189 HU/sec)  
 IEC 60613:2010

X-ray tube assembly cooling is provided by ambient air flow.

X-Ray Tube Assembly  
 Permanent Filtration ..... 0.7mm Al/75 kV IEC 60522/1999

Loading Factors for Leakage Radiation ..... 125 kV, 4.0 mA

Thermal Switch ..... 7A @ 120 Vac or 30 Vdc Max.  
 Normally Closed ..... Opens @ 85°C ±3.9°C (185°F ±7°F)

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




FT 3 Pin - Small diameter 75kV cable with compact taper plug on one end and federal standard terminal on the other.  
 Available in 5 ft. (FT 005) to 100 ft. (FT-100) lengths.

Available through: Parker Medical, Inc.  
 137 New Milford Road East  
 Route 67  
 Bridgewater, CT 06752  
 Phone (203) 350-4304  
 Phone (203) 350-3432

- Notes:
1. Specify straight or 90 degree terminations for both federal and taper connections.
  2. Specify split strain relief for equipment requiring cables to be pulled through small openings.

Stator Ratings and Characteristics  
 Spécificités et Caractéristiques du Stator  
 Statornennleistungen und Merkmale  
 Características y Clarificación de la Bovina

Stator - Wiring Diagram  
 Stator - Schéma de Câblage  
 Stator - Drahtfarbentabelle  
 Bovina - Diagramas

Wire Color Couleurs des Branchements Kabelfarben Cable de Color	Description Description Beschreibung Descripción
1 White Blanc Weiss Blanco	Common Neutre Neutral Común
2 Black Noir Schwarz Negro	Phase Phase Phase Fase
3 Green Vert Grün Verde	Phase Shift Changement de Phase Hilfsphase Cambio de Fase del Estator
4 Red Rouge Rot Rojo	Thermal Switch Switch Thermique Thermoschalter Interruptor Termal
5 Orange Orange Orange Anaranjado	Thermal Switch Switch Thermique Thermoschalter Interruptor Termal
6 Green/Yellow Vert/Jaune Grün/Gelb Verde/Amarillo	Ground Masse Masse Tierra

"R" Stators	Stator "R"	"R" Stator	"R" Bovina	
Black - White	Noir - Blanc	Schwarz - Weiss	Negro - Blanco	16 Ω
Black - Green	Noir - Vert	Schwarz - Grün	Negro - Verde	66 Ω
Green - White	Vert - Blanc	Grün - Weiss	Verde - Blanco	50 Ω
180 Hz Cap	180 Hz Cap	180 Hz Cap	180 Hz Cap	6 μF
60 Hz Cap	60 Hz Cap	60 Hz Cap	60 Hz Cap	30 μF

Stator Drive Frequency Fréquence d'entraînement du stator Statorantrieb Frequenz Frecuencia de la impulsión del estator	RPM
50 Hz	2800 - 3000
60 Hz	3400 - 3600
150 Hz	8500 - 9000
180 Hz	9500 - 10,800

**Stator Power:** 50/60 Hz Only  
 Time to full speed of the anode is a function of the power rating of the "starter" and the weight / diameter of the anode. Time to full speed of 3 inch (71mm) anode series tubes is between 0.5 and 1.0 seconds.

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
### X-Ray Tube Assembly Heating and Cooling Curves

**Note:**  
 Heat inputs into housing include tube power, filament power, and stator power. Heating curves based on no restrictions of natural convection around tube housing assembly.

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Salt Lake City, UT 1-801-972-5000

For a complete listing of our global offices,  
visit [www.vareximaging.com](http://www.vareximaging.com)

Manufactured by Varex Imaging Corporation  
Fabrique par Varex Imaging Corporation  
Hergestellt von Varex Imaging Corporation  
Fabricado por Varex Imaging Corporation

Specifications subject to change without notice.  
Spécifications susceptibles d'être modifiées sans préavis.  
Technische Daten ohne Gewähr.  
Especificaciones sujetas a cambio sin previo aviso.