Note: Document originally drafted in the English language.

**Product Description**

The B-145W housing is designed for Varex Imaging rotating Anode inserts having 71 mm (3 inch) diameter targets.

IEC Classification .....	Class 1
Weight, Approximate:	
Housing & Tube .....	15.74 kg (34.6 lbs)
Mounting .....	Port
Stator Cord .....	3 Feet (0.9m) (Optional)
Cable Receptacles .....	Claymount Mini 75
X-Ray Tube Assembly (Complies to) .....	IEC 60601-2-28

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


Product Description	
Maximum Potential Difference	125 kV
Cathode to Ground	63 kV
Anode to Ground	63 kV
Maximum X-Ray Tube Assembly Heat Content	500 kJ (675 kHU)
Nominal Continuous Input Power	
without Heat Exchanger	140 W (189 HU/sec) IEC 60613:2010
with Heat Exchanger	500 W (675 HU/sec) IEC 60613:2010 (3 liters/min, 20°C - 40°C, 50/50 glycol/water)
Maximum Housing Temperature	78°C
X-Ray Tube Assembly	
Permanent Filtration	0.7mm Al/75 kV IEC 60522/1999
Loading Factors for Leakage Radiation	125 kV, 4.0 mA
Thermal Switch	7A @ 120 Vac or 30 Vdc Max.
Normally Closed	Opens @ 85°C ±3.9°C (185°F ±7°F)
Recirculation Pump	12 VAC, 50/60 Hz, 0.5 A

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

**Stator Power:** 50/60 Hz Only  
 Time to full speed of the anode is a function of the power rating of the "starter" and the weight / diameter of the anode. Time to full speed of 3 inch (71mm) anode series tubes is between 0.5 and 1.0 seconds.

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
### X-Ray Tube Assembly Heating and Cooling Curves

### X-Ray Tube Assembly Heating and Cooling Curves

**Note:**  
 Heat inputs into housing include tube power, filament power, and stator power. Heating curves based on no restrictions of natural convection around tube housing assembly.  
 Cooling performance measured with 50% propylene glycol and water at 3 liters per minute. Use an inhibitor to prevent corrosion of aluminum.

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Salt Lake City, UT 1-801-972-5000

*For a complete listing of our global offices,  
visit [www.vareximaging.com](http://www.vareximaging.com)*

Manufactured by Varex Imaging Corporation  
Fabrique par Varex Imaging Corporation  
Hergestellt von Varex Imaging Corporation  
Fabricado por Varex Imaging Corporation

Specifications subject to change without notice.  
Spécifications susceptibles d'être modifiées sans préavis.  
Technische Daten ohne Gewähr.  
Especificaciones sujetas a cambio sin previo aviso.