Note: Document originally drafted in the English language.

**Product Description**

The B-180H housing is designed for Varex Imaging rotating Anode inserts having 133 mm (5.25 inch) diameter targets.

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Soporte ..... Rotable ó de entrada

### Product Description

|  |                            |
|--|----------------------------|
| Nominal X-ray Tube Voltage .....                             | 150kV                      |
| Maximum Cathode to Ground .....                              | 75 kV                      |
| Maximum Anode to Ground .....                                | 75 kV                      |
| Grid to Cathode  |                            |
| If applicable .....  | -4 kV                      |
| Grid Control Voltages  |                            |
| Typical Bias Voltage for Cutoff at 125 kV .....              | -3600 Vdc                  |
| Grid Voltage for Exposure .....                              | 0 Vdc                      |
| Maximum X-Ray Tube Assembly Heat Content .....               | 1,482 kJ (2,000 kHU)       |
| Nominal Continuous Input Power                               |                            |
| (maximum housing temperature 78°C)                           |                            |
| without HE-100 .....   | 250 W IEC 60613:2010       |
| with HE-100 .....  | 1,130 W IEC 60613:2010     |
| X-ray tube assembly cooling is provided by ambient air flow. |                            |
| X-Ray Tube Assembly  |                            |
| Permanent Filtration .....                                   | 1.0mm Al IEC 60522/1999    |
| Loading Factors for Leakage Radiation .....                  | 150 kV, 7.2 mA             |
| High Voltage Cable Receptacles .....                         | Per IEC 60526              |
| Temperature Limits for Storage and Transport .....           | -20°C to +75°C             |
| Humidity .....   | +10% to +90%               |
| Atmospheric Pressure Range .....                             | 70 kPa to 106 kPa          |
| Thermal Switch .....   | Normally Closed            |
| Open .....   | 76.7°C ±2.8°C (170°F ±5°F) |
|  | 3A @ 120Vac or 30 Vdc max. |
| X-Ray Tube Assembly (Complies to) .....                      | IEC 60601-2-28             |

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






**TRI FOCUS**

**NOTE:**

- REFERENCE DIMENSIONS ARE FOR INFORMATION PURPOSES AND ARE NOT REQUIRED FOR INSPECTION
- APERTURE SIZE OPTIONAL:  
0.937 Sq. in.  
1.250 Sq. in.
- CABLE RECEPTACLE  
2 1/2 -20 UNF -2B THREAD

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
**Stator Power:**  
 Time to full speed of the anode is a function of the power rating of the "starter" and the weight / diameter of the anode. All Varex Imaging stator types are rated for regular speed and high speed starters. Time to full speed of 102 mm (4 inch) anode series tubes is between 1.3 and 2.0 seconds.

Immediately following high speed anode rotation, the rotor speed must be reduced to 4000 r/min or less within 10 seconds using a suitable dynamic braking device.

No more than two high speed starts per minute are permissible. The starting voltage must never exceed 600 volts rms.

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

**X-ray Tube Assembly Heating and Cooling Curve**

**X-ray Tube Assembly Heating and Cooling Curve**

**Caution:**

Housings configured with heat exchangers having long hoses (> 5 ft or 1.5 m) are not configured for operation separate from the associated heat exchanger. In such cases, the tube unit must not be operated unless properly coupled to the heat exchanger. Reference the specific heat exchanger's product data sheet for the associated cooling curves.

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**Note:**

Heat inputs into housing include tube power, filament power, and stator power. Heating curves based on no restrictions of natural convection around tube housing assembly.

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Salt Lake City, UT 1-801-972-5000

For a complete listing of our global offices,  
visit [www.vareximaging.com](http://www.vareximaging.com)

Manufactured by Varex Imaging Corporation  
Fabrique par Varex Imaging Corporation  
Hergestellt von Varex Imaging Corporation  
Fabricado por Varex Imaging Corporation

Specifications subject to change without notice.  
Spécifications susceptibles d'être modifiées sans préavis.  
Technische Daten ohne Gewähr.  
Especificaciones sujetas a cambio sin previo aviso.