Note: Document originally drafted in the English language.

Product Description

The B-199H housing is designed for Varex Imaging rotating anode inserts having 102mm (4 inch) and 120mm (4.75 inch) diameter targets.

IEC Classification Class 1

Weight (approximate)
 Housing & Tube 24 kg (53 lbs)
 With Trunnion Rings 25.5 kg (56 lbs)

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Soporte Rotable ó de entrada

Product Description

Nominal X-ray Tube Voltage 150 kV
 Maximum Cathode to Ground 75 kV
 Maximum Anode to Ground 75 kV

Maximum X-ray Tube Assembly Heat Content 1800 kJ (2,430 kHU)

Maximum Continuous Heat Dissipation
 without air circulator 300 Watts (405 HU/sec)

Maximum Housing Temperature 85°C

X-Ray Tube Assembly
 Permanent Filtration 1.0 mm Al/75 kV IEC 60522

Loading Factors for Leakage Radiation 150kV, 4.0 mA

Temperature Limits for Storage and Transport -20°C to +75°C
 Humidity +10% to +90%
 Atmospheric Pressure Range 70 kPa to 106 kPa

Internal Thermal Switch Normally Closed
 Open 85°C ±3.9°C (185°F ±7°F)
 7A @ 120Vac or 30 Vdc max.

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

Stator Power:

Time to full speed of the anode is a function of the power rating of the “starter” and the weight / diameter of the anode. All Varex Imaging stator types are rated for regular speed and high speed starters. Time to full speed of 102 mm (4 inch) anode series tubes is between 1.3 and 2.0 seconds.

Immediately following high speed anode rotation, the rotor speed must be reduced to 4000 r/min or less within 10 seconds using a suitable dynamic braking device.

No more than two high speed starts per minute are permissible. The starting voltage must never exceed 600 volts rms.

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