Note: Document originally drafted in the English language.

**Product Description**

The DO9 housing is designed for rotating anode inserts.

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Encaje ..... 22.0 kg (48.5 lbs)

Soporte ..... Entrada ó de Rotable

<b>Product Description</b>	
Maximum Potential Difference:	
Cathode to Ground .....	75 kV
Anode to Ground .....	75 kV
Housing Heat Storage Capacity .....	1,200 kJ (1,600 kHU)
Housing Continuous Heat Dissipation	
.....	300 Watts (405 HU/sec)
X-Ray Tube Assembly (includes insert and housing)	
Permanent Filtration .....	0.7mm Al @ 75 kV IEC 60522
Leakage Technique Factors .....	
.....	150 kV, 4.0 mA
Thermal Switch .....	
.....	Normally Closed Rating - 7 A @ 240 VAC
Ambient Air Temperature Limits for Operation .....	
.....	5°C to 40°C
Temperature Limits for Storage and Transportation	
.....	-10°C to 85°C
Humidity .....	10% to 90%
Atmospheric Pressure Range .....	70 kPa to 106 kPa

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


**Stator Power:**  
 Time to full speed of the anode is a function of the power rating of the "starter" and the weight / diameter of the anode. All stator types are rated for regular speed and high speed starters. Time to full speed of 102 mm (4 inch) anode series tubes is between 1.2 and 1.4 seconds.

Immediately following high speed anode rotation, the rotor speed must be reduced to 4000 r/min or less within 10 seconds using a suitable dynamic braking device.

No more than two high speed starts per minute are permissible. The starting voltage must never exceed 600 volts rms.

<b>Starter:</b>	<b>Start</b>	<b>Run</b>
50/60 Hz	220 VAC	50 VAC
150/180 Hz	450 VAC	120 VAC

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