Note: Document originally drafted in the English language.  
Note : Document à l'origine rédigé dans l'anglais.  
Anmerkung: Dokument ursprünglich gezeichnet in der englischen Sprache.  
Nota: Documento elaborado originalmente en la lengua inglesa.

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

**Note:**  
 1. Limits are based on maximum track rating except for the following codes:  
 a - Limited by available heat storage.  
 b - Limited by window heating.  
 c - Limited by filament emission.  
 2. H.S. = Heat Storage  
 kV = Tube Voltage  
 Rating charts reflect maximum tube performance. Tube operation is ultimately limited by system software.

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3 Ø 

|                                                                                                      |                                      |
|------------------------------------------------------------------------------------------------------|--------------------------------------|
| Maximum Peak Voltage .....                                                                           | 140 kV                               |
| Anode to Ground .....                                                                                | 70 kV                                |
| Cathode to Ground .....                                                                              | 70 kV                                |
| Maximum X-ray Tube Assembly Heat Content .....                                                       | 3.8 MJ (5.3 MHU)                     |
| Nominal Continuous Input Power<br>(Includes stator heat) .....                                       | 6.4 kW (9.0 kHU/sec)                 |
|                                                                                                      | IEC 60613:2010                       |
| Maximum Housing Temperature .....                                                                    | 80°C                                 |
| Maximum Heat Exchanger Dissipation .....                                                             | 6.4 kW (9.0 kHU/sec)                 |
| Focal Point Position (Central Ray) Within 1mm<br>(X, Y Direction from the center of radiation port.) |                                      |
| <b>X-Ray Tube Assembly</b>                                                                           |                                      |
| Permanent Filtration .....                                                                           | 0.5 mm Al @ 75 kV IEC 60522/1999     |
| Additional Filtration .....                                                                          | 1.5 mm Al @ 75 kV IEC 60522/1999     |
| Minimum Total Filtration .....                                                                       | 2.0 mm Al @ 75 kV IEC 60522/1999     |
| Loading Factors for Leakage Radiation .....                                                          | 140 kV, 43 mA                        |
| High Voltage Cable Receptacles .....                                                                 | Per IEC 60526                        |
| Ambient Air Temperature Limits for Operation .....                                                   | 5°C to 40°C                          |
| Temperature Limits for Storage and Transport .....                                                   | -20°C to +75°C                       |
| Humidity .....                                                                                       | +10% to +90%                         |
| Atmospheric Pressure Range .....                                                                     | 70 kPa to 106 kPa                    |
| Weight - Housing & Heat Exchanger .....                                                              | 214.5 lbs. ±3.5 lbs.                 |
| IEC Classification .....                                                                             | Class 1                              |
| <b>Safety Devices - Thermal Switches</b>                                                             |                                      |
| External - Normally Closed Contact .....                                                             | Opening at 80°C ±3.9°C               |
| Internal - Normally Closed Contact .....                                                             | Opening at 110°C ±4.4°C              |
| <b>Pressure Switch</b>                                                                               |                                      |
| Normally Open contact .....                                                                          | Contact closes at 3.0 PSID ±1.0 PSID |
| Filament Frequency Limits .....                                                                      | 50 HZ - 40 KHZ                       |
| Power Supply .....                                                                                   | DC                                   |

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





## TUBE HOUSING ASSEMBLY HEATING AND COOLING

**Note:**

- Heat input into housing includes all power sources; tube, filament, stator and circulating pump.
- Heating curves based on no restrictions to air flow through heat exchanger, or natural convection around tube housing assembly.
- Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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

|                         |                |            |
|-------------------------|----------------|------------|
| <b>Stator Type:</b>     |                |            |
| Stator Coil Resistance: |                |            |
| Black to White          | 8.3 Ω ±15%     |            |
| Green to White          | 8.3 Ω ±15%     |            |
| Starter Voltage:        |                |            |
|                         | <u>Start</u>   | <u>Run</u> |
| 105 Hz                  | 480 VAC        | 160 VAC    |
| Starter Current:        |                |            |
|                         | <u>Start</u>   | <u>Run</u> |
| 105 Hz                  | 9 Amps         | 3 Amps     |
| Time to Full Speed:     |                |            |
| 105 Hz                  | 0 - 6300 RPM   | 20 Sec.    |
| X-Ray Tube Assembly:    |                |            |
| GS-5172/B-320H          | IEC 60601-2-28 |            |

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**Note:**  
 Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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