Note: Document originally drafted in the English language.

Product Description

The B-160H housing is designed for Varex Imaging rotating Anode inserts having 108 mm (4.25 inch) diameter targets.

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Soporte Rotable ó de entrada

Product Description

Nominal X-ray Tube Voltage	150 kV
Maximum Cathode to Ground	75 kV
Maximum Anode to Ground	75 kV
Grid to Cathode (If applicable)	-4 kV
Grid Control Voltages	
Typical Bias Voltage for Cutoff at 125 kV	-3700 Vdc Grid
Voltage for Exposure	0 Vdc
Maximum X-ray Tube Assembly Heat Content	1,100 kJ (1,500 kHU)
Nominal Continuous Input Power	
(max. housing temperature 78°C)	
without air circulator	200 Watts (270 HU/sec) IEC 60613:20102010
with air circulator	400 Watts (540 HU/sec) IEC 60613:2010
with HE-100	1000 Watts (1400 HU/sec) IEC 60613:2010
with HE-281	2000 Watts (2800 HU/sec) IEC 60613:2010
with HE-321	2400 Watts (3233 HU/sec) IEC 60613:2010
with HE-581	4000 Watts (5400 HU/sec) IEC 60613:2010
X-ray tube assembly cooling is provided by ambient air flow.	
X-Ray Tube Assembly	
Permanent Filtration	1.0mm Al IEC 60522/1999
High Voltage Cable Receptacles	Per IEC 60526
Loading Factors for Leakage Radiation	
with HE-100	150kV, 6.1mA
with HE-281	150kV, 12.8mA
with HE-321	150kV, 15.4mA
with HE-581	150kV, 15.4mA
Temperature Limits for Storage and Transport	-20°C to +75°C
Humidity	+10% to +90%
Atmospheric Pressure Range	70 kPa to 106 kPa
Thermal Switch	Normally Closed
3 A @ 120 Vac or 30 Vdc max.	
Open	76.7°C ±2.8°C (170°F ±5°F)
X-Ray Tube Assembly (Complies to)	IEC 60601-2-28

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

NOTE:
 1. REFERENCE DIMENSIONS ARE FOR INFORMATION PURPOSES AND ARE NOT REQUIRED FOR INSPECTION
 2. APERTURE SIZE OPTIONAL:
 0.937 Sq. in. 1.250 Sq. in.
 3. CABLE RECEPTACLE
 2 1/2 -20 UNF -2B THREAD

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
Stator Power:

Time to full speed of the anode is a function of the power rating of the "starter" and the weight / diameter of the anode. All Varex Imaging stator types are rated for regular speed and high speed starters. Time to full speed of 102 mm (4 inch) anode series tubes is between 1.3 and 2.0 seconds.

Immediately following high speed anode rotation, the rotor speed must be reduced to 4000 r/min or less within 10 seconds using a suitable dynamic braking device.

No more than two high speed starts per minute are permissible. The starting voltage must never exceed 600 volts rms.

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

Caution:

Housings configured with heat exchangers having long hoses (> 5 ft or 1.5 m) are not configured for operation separate from the associated heat exchanger. In such cases, the tube unit must not be operated unless properly coupled to the heat exchanger. Reference the specific heat exchanger's product data sheet for the associated cooling curves.

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Note:

Heat inputs into housing include tube power, filament power, and stator power. Heating curves based on no restrictions of natural convection around tube housing assembly.

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Salt Lake City, UT 1-801-972-5000

For a complete listing of our global offices,
visit www.vareximaging.com

Manufactured by Varex Imaging Corporation
Fabrique par Varex Imaging Corporation
Hergestellt von Varex Imaging Corporation
Fabricado por Varex Imaging Corporation

Specifications subject to change without notice.
Spécifications susceptibles d'être modifiées sans préavis.
Technische Daten ohne Gewähr.
Especificaciones sujetas a cambio sin previo aviso.