


IHE Statement for Varex Nexus DRF Systems	
Engineering	
Technical Support	
Quality Assurance	

Revision	Date	Originator	Comments
A	2020-3-25	J. Tracy	Initial Version

IHE Integration Statement		Date:	March 2020
Vendor	Product Name		Version
Varex Imaging Corporation	Nexus DRF		3.2 3.5
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below.			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Consistent Presentation of Images (CPI)	Acquisition Modality	None	
Consistent Time (CT)	Time Client	None	
Portable Data for Imaging (PDI)	Portable Media Creator	None	
Scheduled Workflow (SWF)	Acquisition Modality	None	
Radiation Exposure Monitoring (REM)	Acquisition Modality	None	
Varex IHE Statement https://www.vareximaging.com/products/medical/software-solutions/nexus-drf-digital-imaging-system			
Links to Standards Conformance Statements for the Implementation			
Varex DICOM Statement https://www.vareximaging.com/products/medical/software-solutions/nexus-drf-digital-imaging-system			
Links to general information on IHE			
North America: http://www.rsna.org/IHE	Europe http://www.ihe-europe.org	Japan: http://www.jira-net.or.jp/ihe-j	

End of Document