

IEC SYSTEM FOR MUTUAL RECOGNITION OF TEST CERTIFICATES FOR ELECTRICAL EQUIPMENT  
(IECEE) CB SCHEME

**CB TEST CERTIFICATE**

Product Docking Station for x-ray detectors

Name and address of the applicant Claymount Assemblies B.V.  
Anholtseweg, 44  
7091 HB DINXPERLO  
The Netherlands

Name and address of the manufacturer Claymount Assemblies B.V.  
Anholtseweg, 44  
7091 HB DINXPERLO  
The Netherlands

Name and address of the factory

Note: When more than one factory, please report on page 2

Additional information on page 2

Claymount Assemblies B.V.  
Anholtseweg, 44  
7091 HB DINXPERLO  
The Netherlands

Ratings and principal characteristics

Applied part type B, IP20  
Type Harmony 3 (19407-NN, 19408-NN, 19409-NN, 19410-NN) : 12-42 VDC (Fused 1 A@24 V)

Trademark (if any)

VAREX IMAGING

Customer's Testing Facility (CTF) Stage used

Model / Type Ref.

Harmony 3 (19407-NN, 19408-NN, 19409-NN, 19410-NN),  
Harmony 3-ST (20947-NN, 20948-NN & 20949-NN, 20950-NN)  
Harmony 3-NR (19386-NN, 19387-NN, 19388-NN, 19389-NN)

Additional information (if necessary may also be reported on page 2)

Additional information on page 2

where NN - is a 2-digit sequel number which stands for a sequel number to define the exact configuration within the applicable model, e.g. other handle, type of drawer, language etc.  
Re-issue of BE-5986 dated 2015-02-16 due to technical change. New trade mark, update test report to include A1 of IEC 60601-2-54 and IEC 60601-1-2:2014. For details refer to test report.  
The applicant / manufacturer / factory name will change on 01.01.2019 from Claymount Assemblies B.V. to Varex Imaging Nederland B.V.

A sample of the product was tested and found to be in conformity with

IEC 60601-1-3:2008, IEC 60601-1:2005, IEC 60601-1:2005/AMD1:2012, IEC 60601-2-54:2009, IEC 60601-2-54:2009/AMD1:2015

National differences:

CA, CH, US

As shown in the Test Report Ref. No. which forms part of this Certificate

G11N0002, N2290006, N2290007

This CB Test Certificate is issued by the National Certification Body

SGS Belgium N.V. - Division SGS-CEBEC  
Business Riverside Park, Avenue Internationalelaan 55 Build. D, BE-1070 BRUSSEL,  
Belgium